



## Designed for Patients, Physicians and Your office!

Cerene was designed as an in-office global endometrial ablation device. The self-contained system delivers freezing temperatures that effectively ablate endometrial tissue. The small outer-diameter probe, ultra-thin conforming liner and quick treatment time help ensure the procedure is well-tolerated by your patients... without the need for anesthesia or conscious sedation.

### Why consider Cerene?

2.5 minute active ablation time

Self-contained system does not require capital equipment

Minimal anesthesia required – NSAID and/or Paracervical Block


Significant reduction in bleeding episodes and dysmenorrhea symptoms

Post-procedure uterine cavity access and visualization




#### Dysmenorrhea Improvement\*

Percent Reporting Severe or Very Severe Symptoms




#### Median Pain Score\*



#### High Patient Satisfaction\*

with Cerene Treatment




**Indicated use:** Cerene treats heavy menstrual bleeding in premenopausal women due to benign causes for whom child bearing is complete.

# Device Features



Full visualization of uterine cavity in 93% of women 12 mo. post-treatment.\*

\*99% hysteroscopic uterine access.



Hysteroscopic view of uterus at 12-month follow-up visit.

Photo courtesy of Howard Curlin MD and Ted Anderson, MD, PhD, Vanderbilt University Medical Center.



## Cryoablation Benefits

Cryoablation is a well-characterized therapy routinely used to ablate human tissues. In the uterus, the use of cryothermic energy reduces the likelihood of scarring and adhesions while preserving future visualization. It also provides an analgesic effect during treatment.

## Coding and Billing

**CPT 58356** (January 1, 2005)

Private practice reimbursement is dependent on negotiated rates between payers and physicians or physician groups - from general analysis, current payments from private payers and Medicaid patient rates often vary by state. Always check individual payer rates on a case-by-case basis.

Contact us for a demonstration

(510) 338-9301

info@channelmedsystems.com

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Please refer to Instructions For Use for indications, contraindications, warnings and precautions.



**CHANNEL**  
MEDSYSTEMS

*Cool innovation, delivered.*

Channel Medsystems  
2919 7th Street  
Berkeley, CA 94710  
510.338.9301  
channelmedsystems.com  
info@channelmedsystems.com